

# UltraCap® H.D.



## UltraCap® H.D. High Capacity Capsule Filters

Meissner UltraCap® H.D. (Heavy Duty) capsule filters are ready-to-use assemblies that offer high flow and throughput with the convenience and cleanliness of a single-use filter assembly.

Designed for processing of medium to large liquid batches, UltraCap® H.D. assemblies are optimized for continuous and batch processing in biomanufacturing and for final and prefiltration in pharmaceutical, food and beverage, and microelectronics applications. UltraCap® H.D. filters withstand higher operating pressure and are more robust than conventional UltraCap® assemblies.

Meissner UltraCap® H.D. capsule filters are optimized for integration into single-use systems such as Meissner's One-Touch® portfolio.

UltraCap® H.D. assemblies are available with a range of Meissner filter media for liquid, gas, and venting applications. They can be specified with a variety of inlet and outlet connections. An optional gauge port facilitates pressure measurement, while an optional filter stand facilitates fast, easy installation.



### Features and Benefits

- Ruggedized polypropylene assembly withstands higher pressures than conventional high capacity capsule filters and resists damage, ensuring reliability and integrity under demanding conditions
- Encapsulated, integral assembly reduces operator contact with filtered liquids
- Final filtration through prefiltration media options include PVDF, PES, PP and PTFE membranes, as well as PP microfiber, borosilicate glass microfiber, and PP microfiber depth media
- Removal ratings from 0.04 µm to 99 µm
- UltraCap® H.D. filters can be easily configured in series or parallel to maximize design space. 10", 20", 30", 40" and 50" lengths permit fast, easy scale-up
- Extremely low hold-up volume design conserves valuable filtered liquids
- Valved vent port for security and reliability in venting, draining and sampling
- Recessed filter vent/drain on T-style configuration prevents breakage in use
- Single-use filter assembly saves installation, setup, cleaning and cleaning validation costs
- Seamlessly integrates into One-Touch® single-use systems or other single-use portfolios
- Available gamma-irradiated for aseptic applications
- Can be used with UltraSnap™ connectors to configure multiple pre and final capsule filters into a presterilized, ready-to-use assembly

### Materials of Construction

UltraCap® H.D. Housing: Polypropylene (PP)

Filtration Media:

#### Hydrophilic Membranes

SteriLUX® Polyvinylidene fluoride (PVDF)

EverLUX® Polyethersulfone (PES)

STyLUX® Polyethersulfone (PES)

#### Hydrophobic Membranes

Steridyne® Polyvinylidene fluoride (PVDF)

Chemdyne® Polypropylene (PP)

Ultradyn® PTFE

#### Microfiber

ALPHA® Polypropylene (PP)

Vanguard® Polypropylene (PP)

Protec® RF Borosilicate glass (GF)

Protec® RM Borosilicate glass (GF) +  
SteriLUX® PVDF membrane

DeltaMax® Polypropylene (PP) depth

DeltaDepth® Polypropylene (PP) depth

Support Components: Polypropylene (PP)

Sealing Method: Thermal Bonding

### Max. Pressure & Temperature for Liquids

90 psig @ 32°F to 100°F (6,2 bar @ 0°C to 38°C)

55 psig @ 140°F (3,8 bar @ 60°C)

### Max. Pressure & Temperature for Gases

60 psig @ 32°F to 100°F (4,1 bar @ 0°C to 38°C)

35 psig @ 140°F (2,4 bar @ 60°C)

### Connections

Inlet/Outlet: Sanitary flange, hose barb or Flaretek®

Vent Port: Sanitary valve with hose barb

Drain Port: Sanitary valve with hose barb;  
Sanitary plug (T-style option only)

Gauge Port: 3/4" sanitary flange (T-style option only)

### Cartridge Length (Nominal)

10", 20", 30", 40", or 50"

(25 cm, 50 cm, 75 cm, 100 cm, or 125 cm)

### Sterilization

The UltraCap® H.D. assembly must be autoclaved at a minimum of 121°C for 60 minutes with the inlet/outlet down.

UltraCap® H.D. assemblies can be repeatedly autoclaved without loss of integrity.



UltraCap® H.D. assemblies must not be insitu steam sterilized (SIP) as exposure to direct steam flow at 121°C, 15 psig (1 bar) will exceed material design limits and can result in rupture of the plastic housing.

### Mounting

The UltraCap® H.D. assembly can be mounted and supported on suitably braced, rigid, inline pipe connections. A wall mounting bracket and accessory stand are also available. For applications requiring multiple UltraCap® H.D. capsule filters, Meissner's UltraSnap™ filter assembly is recommended. This assembly secures pre and final capsule filters into a single-use filtration system for plug and play use. Contact Meissner for details.



## Configuration Dimensions



## Ordering Information

UltraCap® H.D. Model	Filter Media - Grade		Retention Rating (µm)	Cartridge Length	Body Style	Inlet/Outlet	Vent/Drain Ports
CR2	MF		5	2	T	00	2
CR2 = Standard (non-sterile)	<b>Membrane Media</b>	<b>Grade</b>	<b>Retention Rating (µm)</b>	1 = 10"	T = T-style	00 = 1" sanitary flange	<b>T-Style</b>
	SteriLUX® PVDF	VMH <sup>2</sup> , VTH <sup>1</sup> , VLH <sup>3</sup>	0.1, 0.2, 0.4, 0.6			77 = 3/4" sanitary flange	0 = No vent or drain
	EverLUX® PES	SMH <sup>2</sup> , SLH <sup>3</sup>	0.2, 0.4, 0.6	2 = 20"	N = Inline		1 = No vent; 1/4" sanitary drain plug
GR2 = Gamma-irradiated (sterile)	STyLUX® PES	STW <sup>1</sup> , SLW <sup>3</sup>	0.2	3 = 30"		02 = 1" sanitary flange inlet; 3/8" hose barb outlet	2 = Sanitary vent; 1/4" sanitary drain plug
	STyLUX® PES	SM <sup>2</sup> , SL <sup>3</sup>	0.04, 0.1, 0.2, 0.4, 0.6	4 = 40"		0C = 1" sanitary flange inlet; 1/2" hose barb outlet	3 = Sanitary vent; 3/4" sanitary flange gauge port; 1/4" sanitary drain plug
	STyLUX® PES	ST <sup>1</sup>	0.04, 0.1, 0.2, 0.4	5 = 50"		09 = 1" sanitary flange inlet; 9/16" hose barb outlet	4 = Sanitary vent; no drain
	Steridyne® PVDF	VMV <sup>2</sup> , VTV <sup>1</sup>	0.2			08 = 1" sanitary flange inlet; 3/4" hose barb outlet	5 = Sanitary vent; 3/4" sanitary flange gauge port; no drain
	Gamma-irradiated model not available for media listed below:					0D = 1" sanitary flange inlet; 1" hose barb outlet	6 = No vent or drain; 3/4" sanitary flange gauge port
	Chemdyne® PP	PM <sup>2</sup>	0.04, 0.1, 0.2			22 = 3/8" hose barb	A = No vent; sanitary drain valve
	Chemdyne® PP	PT <sup>1</sup>	0.2			CC = 1/2" hose barb	B = Sanitary vent; sanitary drain valve
	Ultradyn® PTFE	TM	0.05, 0.1, 0.2, 0.4, 1.0, 5.0			99 = 9/16" hose barb	C = Sanitary vent; sanitary drain; 3/4" sanitary flange gauge port
	Ultradyn® PTFE	TA <sup>2</sup> , TT <sup>1</sup>	0.2			88 = 3/4" hose barb	
	<b>Microfiber Media</b>	<b>Grade</b>	<b>Retention Rating (µm)</b>			DD = 1" hose barb	
	ALpHA® PP	MF	0.45, 0.6, 0.8, 1.2, 2.4, 5, 7, 10, 20, 30, 40, 70			AA = 1/2" Flaretek®	
	Vanguard® PP	MN	0.1, 0.2, 0.4, 1, 3, 5, 10, 30, 60, 99 (nominal)			BB = 3/4" Flaretek®	
	Protec® GF	RF*	0.5, 1				<b>Inline</b>
	Protec® GF + PVDF	RM*	0.2, 0.3, 0.5				0 = No vent or drain
	DeltaMax® PP depth	DM	0.5, 1, 3, 5, 10, 20, 40, 70				2 = Two sanitary vent/drain valves
	DeltaDepth® PP depth	DD	0.5, 1, 5, 10, 25, 50 (nominal)				4 = One sanitary vent or drain valve

### Filter Media Grade Descriptions

<sup>1</sup> T-grade (VTH, STW, ST, VTV, PT, TT)

This absolute, microbially rated filter meets full traceability requirements for the pharmaceutical industry. It is 100% integrity tested during manufacture. Each T-grade filter is shipped with a Certificate of Quality stating exact quality control criteria and test performance results. This is a validatable product to meet the stringent requirements of the pharmaceutical industry.

<sup>2</sup> M-grade (VMH, SMH, SM, VMV, PM, TA)

This absolute, microbially rated filter is 100% integrity tested during manufacture. It is suited for critical applications when regulatory documentation requirements are minimal.

<sup>3</sup> L-grade (VLH, SLH, SLW, SL)

This filter is not 100% integrity tested or flushed during manufacture. It is offered as an economical prefilter or final filter when sterility assurance is not required.

SteriLUX®, STyLUX®, EverLUX®, Steridyne®, Chemdyne®, Ultradyn®, ALpHA®, Vanguard®, Protec®, DeltaMax®, DeltaDepth®, UltraCap® and One-Touch® are registered trademarks of Meissner Filtration Products, Inc.

Flaretek is a registered trademark of Entegris Inc.

© 2012, 2007, 2005, 1998 Meissner Filtration Products, Inc. All rights reserved.

CR\_2.0



800.391.9458 • Int'l: +1.805.388.9911 • Camarillo, California  
www.meissner.com